

Fresh Product

Whole Peripheral Blood Vacutainer[®] Tube Collection, Non-Hodgkin's Lymphoma

Catalog#	PB010F-NHL	10 mL
	PB020F-NHL	20 mL
	PB030F-NHL	30 mL
	PB040F-NHL	40 mL
	PB050F-NHL	50 mL
	PB060F-NHL	60 mL
	PB070F-NHL	70 mL
	PB080F-NHL	80 mL
	PB090F-NHL	90 mL
	PB100F-NHL	100 mL

Product Description

Human Whole Peripheral Blood is collected from Non-Hodgkin's lymphoma patients in 10 mL sterile vacutainer(s) containing EDTA.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Blood is obtained using Institutional Review Board (IRB) approved consent forms and protocols.

Sample Collection and Processing

All samples are collected on-site at our Stem Cell Collection Center. Samples are quickly processed in our on-site laboratory to achieve maximum viability and quality.

Format

Fresh whole peripheral blood is shipped at ambient temperature.

Specific containers and media can also be prepared as requested by the customer.

Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

Product Warranty

For our product warranty, please review our Terms and Conditions at stemexpress.com/terms-and-conditions/.

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.