

## Fresh Product

### Peripheral Blood Plasma, Human Immunodeficiency Virus

|          |              |       |
|----------|--------------|-------|
| Catalog# | PBPL005F-HIV | 5 mL  |
|          | PBPL010F-HIV | 10 mL |

## Product Description

Human Peripheral Blood Plasma is collected from human immunodeficiency virus (HIV) patients in 5 mL tubes containing EDTA.

Plasma was obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## Sample Collection and Processing

All samples are collected on-site at our Stem Cell Collection Center. Samples are quickly processed at our on-site laboratory to achieve maximum viability and quality.

## Format

Fresh peripheral blood plasma is shipped at ambient temperature.

## Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

## Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## Product Warranty

For our product warranty, please review our Terms and Conditions at [stemexpress.com/terms-and-conditions/](http://stemexpress.com/terms-and-conditions/).

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.