

## Fresh Product

### Peripheral Blood Serum, Type 1 Diabetes

Catalog#	PBSR005F-DB1	5 mL
	PBSR010F-DB1	10 mL

## Product Description

Human Peripheral Blood Serum from type 1 diabetic patients contains all the components of plasma except for the proteins involved in blood clotting.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Serum was obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## Sample Collection and Processing

All samples are collected on-site at our Stem Cell Collection Center. Samples are quickly processed at our on-site laboratory to achieve maximum viability and quality.

## Format

Fresh peripheral blood serum is shipped at ambient temperature.

## Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

## Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## Product Warranty

For our product warranty, please review our Terms and Conditions at [stemexpress.com/terms-and-conditions/](http://stemexpress.com/terms-and-conditions/).

FOR IN VITRO RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.