

## Fresh Product

### Aneuploidy Maternal Blood Serum Protein Biomarkers Streck<sup>®</sup> Tubes

Catalog#	MBASTR010F-SPB	10 mL
	MBASTR020F-SPB	20 mL
	MBASTR030F-SPB	30 mL
	MBASTR040F-SPB	40 mL
	MBASTR050F-SPB	50 mL
	MBASTR060F-SPB	60 mL

## Product Description

Human Maternal Blood is collected from all types of pregnancies at varying gestations and demographics in 10 mL aliquots.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Maternal blood was obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## Sample Collection and Processing

All samples are collected at our Stem Cell Collection Centers or nearby partner hospitals and clinics. Samples are then quickly processed at our on-site laboratory to achieve maximum viability and quality.

## Format

Fresh maternal blood is shipped at ambient temperature.

## Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

## Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## Product Warranty

For our product warranty, please review our Terms and Conditions at [stemexpress.com/terms-and-conditions/](http://stemexpress.com/terms-and-conditions/).

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.