

## Fresh Product

### Maternal Blood Plasma, Na-Heparin

Catalog#	MBPLNAH005F	5 mL
	MBPLNAH010F	10 mL
	MBPLNAH015F	15 mL
	MBPLNAH020F	20 mL
	MBPLNAH025F	25 mL
	MBPLNAH030F	30 mL

## Product Description

Maternal Blood Plasma is fractionated from vacutainer samples by centrifugation and supplied in 10 mL conical tubes. Harvested plasma contains a reduced platelet count due to the centrifugation process.

Fresh products have a high viability without the detrimental effects of freezing, thawing, and exposure to cryoprotectants.

Plasma was obtained using Institutional Review Board (IRB) approved consent forms and protocols.

## Sample Collection and Processing

All samples are collected at our Stem Cell Collection Centers or nearby partner hospitals and clinics. Samples are then quickly processed at our on-site laboratory to achieve maximum viability and quality.

## Format

Fresh maternal blood plasma is shipped at ambient temperature.

## Storage

Fresh products should be used or processed immediately upon receipt. The warranty only covers items whose specifications are tested at the time they are received.

## Warning

This product contains human tissue or other biological material and MUST be handled at Biosafety Level 2 or higher. All biological products should be treated as potentially infectious or contaminated material, even if infectious disease screening reports are negative. Follow universal precautions and wear appropriate personal protective equipment.

## Product Warranty

For our product warranty, please review our Terms and Conditions at [stemexpress.com/terms-and-conditions/](http://stemexpress.com/terms-and-conditions/).

FOR RESEARCH USE ONLY. NOT APPROVED FOR HUMAN OR VETERINARY USE IN VIVO.